

# COOLGUN



IN 4 MINUTES YOUR FACE WILL CHANGE!

# THE CONCEPT

The Coolifting guns jet sprays the facial tissues using a powerful CO2 flow at a very low temperature and a high pressure, combined with an extremely high concentration of atomized actives. It is the most spectacular and rapid treatment ever seen to challenge the passing of time.

# THE APPLICATION

The combination of the cold CO2 flow at a high pressure and application of atomized anti-aging and lifting ingredients cause:

## Paradoxical vessel-motricity:

Vasoconstriction followed by a vasodilatation of the blood micro-circulation. When suffering an intense cold stimulus, the vessels contract and expand. The removal of localised waste products is expedited.

## Thermal Shock:

Stimulation of the underlying tissues. Those tissues are normally at a temperature of some 36°C (97°F). The Coolifting extreme and sudden cold causes a deep and instantaneous reaction in the dermis.

## Cryophoresis and Barophoresis:

The combination of cold accompanied with the high pressure applied on the dermis, causes an intense penetration of the Coolifting actives and an expanding increase of their effects.

## Cryo-carboxytherapy:

The CO2 flow at a very low temperature and at a high pressure activates the blood micro-circulation; it induces collagen and elastin formation and increases the tissues' oxygenation.



# ABOUT COOL LIFTING

It is a revolutionary registered system that allows performing a shocking facial treatment through a safe, fast and efficient method.

## How does it work?

The Coolifting gun jet sprays the facial tissues using a powerful CO2 flow, at a very low temperature and a very high pressure, combined with a high concentration of atomized anti-aging and lifting ingredients.

## How many sessions do I need?

One session (4 Mins): Visibly Smoother, nourished and more radiant skin. We recommend 6 treatments in close succession to roll back up to 10 years from your skin.

## What will I feel during and after the session?


The procedure is painless. The application, pressure and temperature are technically controlled and doesn't involve any discomfort. No special care is required following the session.




BEFORE



AFTER



BEFORE



AFTER

Visible effects in just one session.

Intense hydration-collagen stimulation-dermis regeneration-improvement of skin texture and quality.



0800 3100 421 | [info@cryosthetics.co.uk](mailto:info@cryosthetics.co.uk)  
[www.cryosthetics.co.uk](http://www.cryosthetics.co.uk)